



SANJMED MEDICAL DISTR. CC.

BP MONITORING



Specifications:

- Economy type blood pressure monitor
- Shock-proof chrome plated gauge with protective ring
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Red zone on gauge between 80mmHg and 160 mmHg
- Navy blue adult size cotton cuff
- 2-tube adult size cuff and bladder
- Standard latex bulb and inflation valve with spring
- Vinyl zipper case



Specifications:

- Square shaped 14cm diameter dial
- 360° Manometer swivel for easy-reading
- Left and right adjustable
- Aluminum stand and tube
- 5 wheel base with brake
- Adjustable height (range of 850mm to 1250mm)
- Metal basket at rear of unit
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Red zone on gauge between 80mmHg and 160 mmHg
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder
- PVC coil tubing
- Standard latex bulb and inflation valve with spring



SANJMED MEDICAL DISTR. CC.



Specifications:

- Square shaped 145 mm diameter dial
- Front and rear adjustable for easy reading

Specifications:

- Square shaped 140mm diameter dial
- 90° Manometer swivel for easy-reading
- Left and right adjustable
- Metal basket behind unit
- Aneroid gauge calibrated from 0mmHg to 300mmHg
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder
- PVC coil tubing
- Standard latex bulb and inflation valve with spring

lg
g

Specifications:

- Glass tube (3mm) with graduation scale
- Mercurial gauge calibrated from 0mmHg to 300mmHg
- Heavy chrome-plated stainless steel case
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder with D-ring
- Standard latex bulb and inflation valve with spring



Specifications:

- Glass tube (3mm) with graduation scale
- Mercurial gauge calibrated from 0mmHg to 300mmHg
- Red zone on gauge between 80mmHg and 160 mmHg
- Navy blue washable adult size nylon cuff with D-ring
- 2-tube adult size cuff and bladder
- Standard latex bulb and inflation valve with spring
- PVC coil tubing
- Aluminum stand and tube
- 5 wheel base with brake
- Adjustable height (range of 850mm to 1250mm)

